





Medi-3200

Medical Imaging Systems



Professional Mammography Film Digitizer



Introduced by Microtek (ISO13485 certified), the professional design of Medi-3200 can easily transfer mammography film into digital format, simply using customized templates and the One-Pass-Scan feature together, the Medi-3200 can easily carry out the auto-crop function to scan multi-films at a time. In addition, it captures details in bright and dark areas at a high optical resolution up to 3200 dpi, 4.2 Dmax, and 16-bit grayscale. Therefore, it provides the medical professionals a smart and easy way to digitize and observe the mammography films.

Furthermore, the digitized mammography films can be easily read and utilized in the optional DICOM gateway software to be integrated with HIS/PACS system.

Advantages

- The most economical mammography film digitizer
- One-Pass-Scan to auto-crop multi-films scan at a time
- Provides high optical resolution up to 3200 dpi for better image detail
- Offers dynamic range up to 4.2 Dmax
- Covers various sizes of x-ray films up to 12"x16"
- Supports telemedicine with an optional DICOM software
- Minimizes the space and cost of the film storage
- Offers better quality, excellent performance, and secures records