




# MiDICOM



Medical Imaging Systems



## DICOM Converter Solution

Microtek's MiDICOM provides a reliable converter solution for bringing texts and images together to create DICOM-format medical reports. By applying with MiDICOM, varies e-files formats, such as Microsoft Word, Excel, Powerpoint, Adobe PDF, and JPEG, can be transformed into DICOM format in seconds. These DICOM-covered files can be directly previewed and uploaded to PACS/RIS systems, quickly and efficiently managed and shared among the medical professionals.

In addition, working MiDICOM together with Microtek's scanners, users can tremendously reduce paper and chemical waste associated with standard physical modalities, which means no more fuzzies to deal with physicians' handwritings, diagnostic reports or any hard-copy files. All files will be electronicalized into DICOM format quickly.

## Features

- Downloads data directly from the worklist server
- Allows to create patient information manually
- Converts varies file formats easily into DICOM format
- Uploads data to PACS/RIS server